

Prepare cells as per usual procedure but when required to aliquot 50 μ l of cells into tubes, aliquot into wells of 96 well U-bottom or V-bottom trays instead (NOT flat-bottomed as cells won't pellet during centrifugation – V-bottom trays are best)

Staining procedure is exactly the same as when performed in tubes - the same ratio of antibody : cells will be required, however wash volumes can be reduced to 200 μ l / wash and still only one wash required between antibody incubation steps.

Wash technique:

- Add 200 μ l of PBA to each well
- Centrifuge – 350 x *g* for 1 minute
- Gently flick out buffer with minimum force required to prevent losing cell pellet
- With plate right-side up, gently tap plate on sides to resuspend cell pellets in buffer remaining in wells

UNCONTROLLED DOCUMENT