

## Flow Cytometry (murine): Direct Protocol

Please note – this protocol has been optimised for certain conditions and is meant as a guide only, you may need to alter some parameters to suit your own experiment

1. Wash cells in PBS, and resuspend to approx  $4 \times 10^6$  viable cells/ml in PBS/BSA/azide
2. Add 5  $\mu$ l of murine gamma-globulin (at 10 mg/ml) per  $10^6$  cells, to block Fc receptors.
3. Incubate at room temperature for at least 30 minutes, then place on ice. There is no need to wash out the gamma-globulin.
4. Aliquot 50  $\mu$ l of Fc-blocked cells into ice-cold FACS tube (Falcon 2008)
5. Add 10  $\mu$ l of each primary antibody and mix. For all primary antibodies the optimal concentration should be determined by titration. Each antibody must be directly labelled with a fluorochrome or biotin, and all antibodies that are to be used together in one tube must be labelled with *different* conjugates.
6. Incubate on ice 30 minutes (*during incubation, set centrifuge to cool to 4°C*)
7. Add 3 ml PBS/BSA/azide, spin down, tip off supernatant
8. If one antibody is biotinylated:
  - dilute streptavidin-conjugate to a suitable concentration in PBS/BSA/azide
  - add 50  $\mu$ l to cells.
  - incubate on ice for 30 minutes.
9. Wash with 3 ml ice-cold PBS/BSA/azide.
10. Wash all tubes with 3ml PBS/azide (NO BSA).
11. Resuspend in 200 – 500  $\mu$ l paraformaldehyde (depending on number of cells) to fix
12. Store at 4°C in the dark until analysis (up to 1 week).

### Suggested controls\*:

1. Cells only
2. Isotype-matched negative control antibodies for each fluorochrome, or FMO controls
3. Compensation controls (ie, a brightly-staining antibody by itself) for each fluorochrome

**\*NOTE:** See Recommended Controls document for comprehensive details