

Labelling Human Cells for Flow Cytometry: Indirect Protocol (1-3 colour)

Please note – this protocol has been optimised for certain conditions and is meant as a guide only, you may need to alter some parameters to suit your own experiment

1. Wash cells at least once and resuspend to approx 4×10^6 viable cells/ml in staining buffer.
2. Add purified human IgG to block Fc receptors. Use 5 μ l of 10 mg/ml stock per million cells. Incubate at room temperature for 30 minutes
3. Aliquot 50 μ l of Fc-blocked cell suspension into FACS tubes (Falcon 2008)
4. Add 10 μ l of unlabelled primary antibody and mix. For all primary antibodies the optimal concentration should be determined by titration. If also using a directly-labelled antibody/s (ie multicolour staining), **DO NOT** add at this step.
5. Incubate at room temperature for 30 minutes (*during incubation, set centrifuge to cool to 4°C and dilute secondary antibody - see Step 7. Also dilute streptavidin conjugate, if required - see Step 10*)
6. Add 3 ml ice-cold staining buffer, spin down, tip off supernatant (*allow to drain briefly after tipping off, and be consistent between tubes, as different amounts of liquid left in bottom of tube may affect subsequent steps*)
7. Dilute secondary antibody (eg anti-IgG (or anti-IgM) conjugated to a fluorochrome or biotin) to appropriate dilution in staining buffer, and add 1% human IgG. Allow to pre-adsorb for at least 20 minutes at room temperature.
8. Add 50 μ l of pre-diluted secondary antibody to tubes and vortex for 3 seconds. Incubate on ice 30 minutes
9. Add 3 ml ice-cold staining buffer, spin down, tip off supernatant, vortex for 3 seconds
10. If doing multicolour labelling:
 - add 10 μ l of purified IgG corresponding to the host of the primary antibody at 2 mg/ml per tube (this is CRITICAL to block free binding sites on the secondary)
 - incubate for 20 minutes at room temperature
 - add directly labelled primary antibody/s (10 μ l at 50 μ g/ml)
 - if secondary antibody is biotinylated instead of directly labelled, then also add 50 μ l of streptavidin-fluorochrome conjugate (diluted appropriately, with 1% hulgG if appropriate)
 - incubate at room temperature for 30 minutes
 - add 3 ml ice-cold staining buffer, spin down, tip off supernatant, vortex for 3 seconds
11. Add 3 ml ice-cold PBS/azide (NO serum) to ALL tubes. Spin down, tip off supernatant
Resuspend in 200 – 500 μ l 1% paraformaldehyde and store at 4°C in the dark until analysis

Suggested controls*:

1. Cells only
2. Isotype-matched negative control antibodies for each fluorochrome, or FMO controls
3. Compensation controls (ie, a brightly-staining antibody by itself) for each fluorochrome

***NOTE:** See Recommended Controls document for comprehensive details